

ReadFR - SNP Quality Control Program

User Manual and Technical Guide

Program Version: 1.0 - GenomeQC SNP Quality Control Pipeline

Release Date: February 13, 2026

Author: Dr. DEUKMIN LEE

Organization: Hankyong National University

Department: Department of Animal Science

Email: dhlee@hknu.ac.kr

Table of Contents

1. Overview
 2. System Requirements
 3. Installation
 4. Quick Start Guide
 5. Parameter File Configuration
 6. Input Data Formats
 7. Quality Control Criteria
 8. Output Formats
 9. Advanced Usage
 10. Troubleshooting
 11. Case Studies
 12. References
-

1. Overview

What is ReadFR?

ReadFR is a high-performance Fortran program designed for comprehensive quality control (QC) of SNP genotype data from Illumina FinalReport files. It processes large-scale genomic datasets, applies multiple QC criteria, and generates standardized output formats compatible with GBLUP and other genomic analysis tools.

Key Capabilities

- **Scale:** Processes 10,000+ SNPs and 1,000+ animals simultaneously
- **Speed:** $O(1)$ hash table-based animal lookup
- **Flexibility:** Configurable QC thresholds via parameter file
- **Output:** BLUPF90-compatible GENO format
- **Reliability:** Comprehensive error checking and reporting

Main Features

Feature	Description
Hash Table Technology	O(1) animal lookup time
Memory Efficiency	Dynamic memory allocation
QC Filters	GC Score, R-Intensity, GT Score, Cluster Separation, Call Rate
Multi-parameter Support	Multiple SNP chips and versions
Flexible Configuration	Parameter file-based settings
Detailed Reporting	QC statistics and filtering results

2. System Requirements

Hardware Requirements

Component	Minimum	Recommended
CPU	2 cores	4+ cores
RAM	4 GB	8 GB+
Storage	100 MB	500 MB+
Network	Optional	For data transfer

Software Requirements

Software	Minimum Version	Recommended
OS	Linux (CentOS 7+)	Ubuntu 20.04 LTS+
Compiler	gfortran 4.8	gfortran 9.0+
Tools	make, ar	make, ar, gdb

Data Requirements

- **FinalReport File:** Illumina SNP genotyping output
- **Pedigree File:** Animal information and relationships
- **MAP File:** SNP physical positions and marker information
- **DATA File** (optional): Additional animal information

3. Installation

Quick Installation (System-wide)

```
cd DKBLUPF90
sudo ./install.sh
```

User Directory Installation (No root required)

```
cd DKBLUPF90
PREFIX=$HOME/.local ./install.sh
```

```
export PATH="$HOME/.local/bin:$PATH"
export LD_LIBRARY_PATH="$HOME/.local/lib:$LD_LIBRARY_PATH"
```

Verification

```
ReadFR --help
which ReadFR
```

See **INSTALL.md** for detailed installation instructions.

4. Quick Start Guide

Step 1: Prepare Parameter File

Create a parameter file (named `parameter`):

```
# SNP File information
SNPFILE: PorcineSNP60_SJ_51th_595sp_DY_FinalReport.txt
HEADER: 10
DELIM: TAB
NO_VARIABLES: 11
2 ANIMAL_ARN
5 SNP_ID
9 CHR
10 POS
11 ALLELE1_AB
12 ALLELE2_AB
25 R_INTENSITY 0.4 2.0
27 GC_SCORE 0.65
30 GT_SCORE 0.50
31 CLUSTER_SEP 0.30
99 CALL_RATE 0.70
```

```
# MAP File
MAPFILE: MAP_K.txt
HEADER: 1
DELIM: TAB
NO_VARIABLES: 5
2 SNP_ID
3 CHR
4 POS
5 ARRAY_ALL
6 ARRAY_CHR
```

```
# PED File
PEDFILE: PED_Total.txt
HEADER: 0
DELIM: TAB
```

```
NO_VARIABLES: 7
```

```
1 BREED
```

```
2 ID
```

```
3 ARN
```

```
4 SIRE
```

```
5 DAM
```

```
6 SEX
```

```
7 BDATE
```

```
# Output prefix
```

```
OUTPUTPREFIX: GENO_QC
```

Step 2: Prepare Input Files

1. **FinalReport File:** Export from GenomeStudio
2. **MAP File:** SNP information
3. **PED File:** Animal pedigree information

Step 3: Run Program

ReadFR parameter

Step 4: Check Output

```
ls -lh GENO_QC_*.geno
```

```
head -2 GENO_QC_YYYYMMDD_00.geno
```

5. Parameter File Configuration

Basic Structure

```
# Comments start with #
```

```
KEYWORD: value1 value2 ...
```

File Keywords

SNPFILE Illumina FinalReport file path

```
SNPFILE: path/to/FinalReport.txt
```

PEDFILE Pedigree file path

```
PEDFILE: path/to/pedigree.txt
```

MAPFILE SNP map information file path

```
MAPFILE: path/to/snp_map.txt
```

DATAFILE (Optional) Additional animal information file

DATAFILE: path/to/additional_data.txt

OUTPUTPREFIX Prefix for output files

OUTPUTPREFIX: GENO_QC

Output: GENO_QC_YYYYMMDD_00.geno, GENO_QC_YYYYMMDD_01.geno, ...

File Format Keywords

HEADER Number of header lines to skip

HEADER: 10 # Skip first 10 lines

DELIM Delimiter type: TAB, SPACE, or comma

DELIM: TAB # Tab-delimited

DELIM: SPACE # Space-delimited

NO_VARIABLES Number of variables/columns to use

NO_VARIABLES: 11 # Use first 11 columns

Field Mapping

Each field is defined by: column_number FIELD_NAME [threshold1 threshold2]

QC Thresholds Some fields accept additional threshold values:

25 R_INTENSITY 0.4 2.0 # min=0.4, max=2.0

27 GC_SCORE 0.65 # min=0.65

30 GT_SCORE 0.50 # min=0.50

31 CLUSTER_SEP 0.30 # min=0.30

99 CALL_RATE 0.70 # min=0.70

Complete Example

```
# =====  
# SNP File Configuration  
# =====  
SNPFILE: GenomeStudio_Export.txt  
HEADER: 10  
DELIM: TAB  
NO_VARIABLES: 11  
2 ANIMAL_ARN  
5 SNP_ID  
9 CHR  
10 POS  
11 ALLELE1_AB  
12 ALLELE2_AB
```

```

25 R_INTENSITY 0.4 2.0
27 GC_SCORE 0.65
30 GT_SCORE 0.50
31 CLUSTER_SEP 0.30
99 CALL_RATE 0.70

# =====
# MAP File Configuration
# =====
MAPFILE: SNP_positions.txt
HEADER: 1
DELIM: TAB
NO_VARIABLES: 5
2 SNP_ID
3 CHR
4 POS
5 ARRAY_ALL
6 ARRAY_CHR

# =====
# Pedigree File Configuration
# =====
PEDFILE: Animal_pedigree.txt
HEADER: 0
DELIM: TAB
NO_VARIABLES: 7
1 BREED
2 ID
3 ARN
4 SIRE
5 DAM
6 SEX
7 BDATE

# =====
# Output Configuration
# =====
OUTPUTPREFIX: Analysis_Result

```

6. Input Data Formats

FinalReport File Format

Illumina GenomeStudio export (tab-delimited):

```

Header Line 1
Header Line 2

```

...

Header Line 10 (usually contains column names)

Animal_ARN	X-Coordinate	Y-Coordinate	...	SNP_ID	...	R-Intensity	GC_Score	GT_Score	Clus
ARN001	12345	67890	...	snp_1	...	1.85	0.78	0.95	0.4
ARN002	11234	68901	...	snp_1	...	1.92	0.81	0.92	0.4

MAP File Format

SNP position information (tab-delimited):

SNP_ID	CHR	POS	ARRAY_ALL	ARRAY_CHR
snp_1	1	12345678	1	1
snp_2	1	23456789	2	2
snp_3	2	34567890	3	1

PED File Format

Animal pedigree information (tab or space-delimited):

BREED	ID	ARN	SIRE	DAM	SEX	BDATE
Duroc	PK001	ARN001	0	0	2	20200101
Duroc	PK002	ARN002	ARN001	0	1	20210315

7. Quality Control Criteria

GC_SCORE Filter

Purpose: Illumina genotyping quality metric

Valid Range: 0.0 - 1.0

Recommended Threshold: 0.65

R_INTENSITY Filter

Purpose: Overall signal intensity

Valid Range: 0.0 - 3.0

Recommended Range: 0.4 - 2.0

Note: Requires TWO values: min and max

GT_SCORE Filter

Purpose: Genotype clustering quality

Valid Range: 0.0 - 1.0

Recommended Threshold: 0.50

CLUSTER_SEP Filter

Purpose: Cluster separation quality

Valid Range: 0.0 - 1.0

Recommended Threshold: 0.30

CALL_RATE Filter

Purpose: Animal-level call rate

Valid Range: 0.0 - 1.0

Recommended Threshold: 0.70

Default Thresholds

Criterion	Default Value	Recommendation
GC_SCORE	0.65	0.65-0.75
R_INTENSITY min	0.4	0.3-0.5
R_INTENSITY max	2.0	1.8-2.2
GT_SCORE	0.50	0.50-0.60
CLUSTER_SEP	0.30	0.25-0.35
CALL_RATE	0.70	0.90-0.95

8. Output Formats

Output File Naming

[PREFIX]_[YYYYMMDD]_[SequenceNumber].geno

Example: GENO_QC_20260213_00.geno

GENO File Format

Header line:

Animal_ID BREED SIRE DAM SEX BDate LOC GENO(1-76756)

Data lines:

PK001 Duroc 0 0 2 20200101 DY 0 1 1 0 0 1 2 -1 ...
PK002 Duroc ARN001 0 1 20210315 DY 2 0 0 1 1 0 1 9 ...

GENO File Specifications

- **Genotype Coding:** 0 (homozygous 1/1), 1 (heterozygous 1/2), 2 (homozygous 2/2), 9 (missing)
 - **Column Order:** Fixed (as shown in header)
 - **Format:** Space-delimited
 - **Encoding:** Plain text, gzip-compressed available
-

9. Advanced Usage

Case-Insensitive Parameter Handling

The program accepts parameter files with flexible casing:


```
# All accepted formats:
snpfile: data.txt
SNPFILE: data.txt
SnpFile: data.txt
SNPFile: data.txt
```

Field Name Flexibility

Field names support multiple formats:

```
# All equivalent:
25 R_INTENSITY 0.4 2.0
25 R-INTENSITY 0.4 2.0
25 r_intensity 0.4 2.0
25 r-intensity 0.4 2.0
```

Large Dataset Processing

For datasets with >100,000 animals:

```
# Increase stack memory
ulimit -s unlimited

# Run with monitoring
time ReadFR parameter

# Check memory usage
top -p $(pgrep -f "ReadFR parameter")
```

Parallel Batch Processing

Process multiple files sequentially:

```
for file in *.txt; do
    echo "Processing $file..."
    sed "s|SNPFILE:.*|SNPFILE: $file|" parameter.template > parameter_${file}
    ReadFR parameter_${file}
    mkdir -p results_${file}
    mv GENO_QC_*.geno results_${file}/
done
```

10. Troubleshooting

Issue: “Input parameter file” error

Cause: No parameter file provided

Solution:

```
ReadFR parameter
# Not: ReadFR
```

Issue: “cannot open file” for input data

Cause: File path error or file not found

Solution:

```
# Check file exists  
ls -l /path/to/file  
  
# Use absolute paths in parameter file  
# Relative paths are OK from same directory
```

Issue: Segmentation fault during execution

Cause: Memory issue or corrupted data

Solution:

```
# Check available memory  
free -h  
  
# Reduce dataset size for testing  
# Check file integrity  
file yourfile.txt
```

Issue: Unexpected filtering results

Cause: Incorrect threshold values or data format

Solution:

```
# Verify thresholds in parameter file  
# Check data value ranges in source files  
# Test with default thresholds first
```

11. Case Studies

Case Study 1: Large-Scale Commercial Genotyping

Scenario: 1000 pigs, 60K SNPs

Processing Time: ~5 minutes

Output Size: 15 MB

Parameter Configuration:

```
SNPFILE: Commercial_GenomeStudio_Export.txt  
SNPFILE: Commercial_Pedigree.txt  
MAPFILE: SNP60k_map.txt  
OUTPUTPREFIX: Commercial_QC
```

```
# Strict QC for downstream analysis  
GC_SCORE: 0.75  
R_INTENSITY: 0.5 1.9  
CALL_RATE: 0.95
```

Case Study 2: Research Project with Multiple Chips

Scenario: Mixed SNP50K and SNP60K

Processing Approach: Separate parameter files per chip

SNP50K_parameter:

SNPFILE: Chip_50K_FinalReport.txt

MAPFILE: SNP50k_map.txt

OUTPUTPREFIX: Analysis_50K

SNP60K_parameter:

SNPFILE: Chip_60K_FinalReport.txt

MAPFILE: SNP60k_map.txt

OUTPUTPREFIX: Analysis_60K

12. References

Related Publications

- Lee, D. (2025). “Hash table-based genome data processing for large-scale genomic selection.” Journal of Genomics.
- GenomeStudio User Guide. Illumina, Inc.
- BLUPF90 Documentation. University of Georgia.

External Resources

- Illumina Genotyping
- BLUPF90 Suite
- Fortran Documentation

Support

For technical support: - Email: dhlee@hknu.ac.kr - Institution: Hankyong National University - Department: Department of Animal Science

Appendix: Version History

Version	Date	Changes
1.0	Feb 13, 2026	Initial release

License

DKBLUPF90 and ReadFR are released under the MIT License.

Document Information - Title: ReadFR SNP Quality Control Program - User Manual - Version: 1.0 - Date: February 13, 2026 - Author: Dr. DEUKMIN LEE - Organization: Hankyong National University

End of Document